

On-Call[®] Express

Blood Glucose Monitoring System

Reliable

- Clinically proven accuracy
- Auto calibration circuitry in every test strip
- Audible fill detection and redosing to avoid wasted strips

Fast

- 4 second testing

Less Pain

- Tiny 0.4 µL blood sample
- Alternate site testing (Palm, Forearm)

Easy

- No coding required
- User friendly design
- Large display with bold numbers
- Convenient test strip ejector

Advanced Features

- Stores 300 tests with date/time
- 7, 14, 30, 60 and 90-day averages
- Helpful meal markers
- Hypo and hyper alarms
- Ketone test alert
- 5 daily testing reminders



Feature	Specification
Measurement Range	20 - 600mg/dl
Result Calibration	Plasma-equivalent
Minimum Sample Size	0.4 µL
Test Time	4 seconds
Battery Life	1,000 tests
Glucose Units of Measure	The meter is pre-set to milligrams per deciliter (mg/dl)
Memory	Up to 300 records with time and date
Automatic shutoff	2 minutes after last action or 5 sec after strip removal
Meter Size	85.0 mm x 49.0 mm x 16.4 mm
Display Size	35 mm x 32 mm
Operating Temperature	41-113°F (5-45°C)
Operating Relative Humidity	10 - 90% (non-condensing)
Hematocrit Range	30-55%
Data Port	USB

24/7 Support: 1-800-838-9502

www.oncallmeters.com



©2015 ACON Laboratories, Inc. San Diego, CA 92121
PN 2150009302 6/15