



H. pylori Ag Cassette FAQ Frequently Asked Questions (FAQ)

Q: How does the test work?

A: The *H. pylori* Antigen Rapid Test Cassette (Feces) is a qualitative, lateral flow immunoassay for the detection of *H. pylori* antigens in human feces specimens. In this test, the membrane is pre-coated with anti-*H. pylori* antibodies on the test line region of the test cassette. The mixture migrates upward by capillary action to react with anti-*H. pylori* antibodies on the membrane and generate a colored line. The presence of this colored line in the test line region indicates a positive result, while its absence indicates a negative result. To serve as a procedural control, a colored line will always appear in the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

Q: What can this test help diagnose?

A: This test will help diagnose *H. pylori* infection. This infection may cause:

- Ulcers
- Inflammation of the Stomach Lining
- Stomach Cancer

Q: What is the difference between an antibody test and an antigen test for *H. Pylori*?

A: The antibody test for *H. pylori* is performed on serum, plasma, or whole blood. This test would detect antibodies from a current or previous infection. The *H. pylori* antigen test is performed on feces and would only detect a current infection.

Q: How accurate is the test?

A: The *H. pylori* Antigen Rapid Test Cassette (Feces) has been compared with Elisa methods, demonstrating an overall accuracy of 99%. It features a relative sensitivity of 99% and relative specificity of 98.9%.

Q: Are there any patient restrictions?

A: Following certain antibiotic treatments, the concentration of *H. pylori* antigens may decrease to the concentration below the minimum detection level of the test. Therefore, diagnosis should be made with caution during antibiotic treatment.

Q: How do I know that the test was run properly?

A: The appearance of a colored line at the Control Line region tells you that you followed the test procedure properly and the proper amount of sample was used.



aconlabs.com

ACON Laboratories, Inc.
10125 Mesa Rim Road, San Diego, CA 92121, U.S.A.
Tel: 1.858.875.8000
Fax: 1.858.200.0729
Email: info@aconlabs.com