



MISSION®
**CHOLESTEROL
METERS**

Excellence

in Diagnostic Products

About ACON



San Diego, California USA

ACON Laboratories, a global, award winning manufacturer with its headquarters in San Diego, has proudly built its reputation in the medical device industry as a strategic partner and OEM manufacturer for some of the largest companies worldwide. At the core of the organization resides the adage, “Your Trusted Partner.” The company’s goal is to make a positive health impact on the millions of people around the world that use its products.

ACON Laboratories found initial success by revolutionizing the lateral flow rapid test industry, becoming a market leader in the in vitro diagnostics industry. Since then, it has increased its point-of-care product portfolio to include Diabetes care, Urinalysis, Cholesterol, Coagulation, Hemoglobin and Chemistry tests.

Its manufacturing campus in Hangzhou, China is recognized as a state-of-the-art facility for its manufacturing operation and production of industry-leading products. The facility is ISO 13485:2003 certified and was the first rapid diagnostics manufacturing facility in China to be FDA registered. For three consecutive years, it was recognized by Forbes China Top 100. Today, ACON Laboratories supplies over 1 billion products annually across 140+ countries.

Certificates and Registrations

- ▶ FDA registered and inspected CE Mark
- ▶ cGMP/QSR compliant (Class 100,000 clean room for manufacturing)
- ▶ ISO 13485
- ▶ ISO 15197
- ▶ IVDD and MDD certified

ACON Timeline

1996

ACON was founded in Bethlehem, PA, USA. ACON operated in a 27,000 sq. ft. manufacturing facility.

2001

Increased its manufacturing facility to 150,000 sq. ft. with more than 1,500 employees.

2009

Expanded product lines. ACON sold the LF rapid diagnostic business in Asia, the Middle East, Africa and Latin America. ACON continued to expand its business in Diabetes Care, Clinical Chemistry, and Immunoassay.

2016

Launched the PT/INR coagulation monitoring system in Europe and Latin America.

Your **Trusted** Partner

ACON works strategically with its sister companies to produce high quality medical products that meet and exceed the customers exacting quality expectations.

We are committed to continuously monitoring our quality systems to ensure their effectiveness and compliance to all applicable regulatory and quality requirements. Exceeding customer expectations and maintaining customer loyalty are our top priorities.

1999

Moved to San Diego, California, USA

2006

Launched new product lines. ACON sold the lateral flow (LF) rapid diagnostic business in the US, Europe, Canada, Israel, Japan, Australia, and New Zealand. ACON launched new product lines including Diabetes Care, Clinical Chemistry, and Immunoassay.

2014

Expanded and moved to a larger facility for manufacture, R&D and business development. Our new 70,000 sq. m (750,000 sq. ft.) facility includes state of the art manufacturing equipment to supply the growing demand of the global diagnostic market. The new facility will also allow our R&D and business development teams to expand, enabling ACON to develop more innovative products to meet our customer's needs worldwide.

2018


Launched the On Call[®] Care Alliance platform allowing healthcare providers and care teams to track, manage, and engage patients remotely.



The **Mission Cholesterol Meters** and test strips are known throughout the world for their quality, ease of use and outstanding performance.

When using Mission® Cholesterol meters, one can expect only accurate results. Ergonomically designed with the end user in mind, the Mission Cholesterol meters from ACON exudes confidence, accuracy, durability and a peace of mind that you would expect from a leading company in diagnostic instrument.

The outstanding growth in Mission Cholesterol meters and strips can only be contributed to the confidence and trust our customers have in these products.



Lipids

Mission Ultra Cholesterol

SPECIFICATIONS & FEATURES

- 4 μ L Capillary whole blood
- Total CHOL results in as little as 15 seconds
- Biosensor technology for direct sampling application
- Memory for 200 results
- Operated in battery or standard AC adaptor mode



Mission Cholesterol Meter

SPECIFICATIONS & FEATURES

- 10 μ L (individual test), 35 μ L (3-in-1 lipid panel test)
- Results in \leq 2 minutes
- Capillary or venous whole blood, serum or plasma
- Accuracy comparable to laboratory analyzers
- Reflectance Photometry



MISSION CHOLESTEROL 3 IN 1 LIPID TEST 3 IN 1

- Tests for the complete lipid panel
 - Total Cholesterol (CHOL)
 - High-Density Lipoprotein (HDL)
 - Triglycerides (TRIG)
 - Calculated Low-Density Lipoprotein, CHOL/HDL ratio
 - Cardiac Risk Assessment (CHD)
- 18-months extended shelf life
- Room temperature storage
- Individually pouched test devices preserve stability



MISSION CHOLESTEROL METER FEATURES

- Large, easy-to-read display
- Fast, reliable lipid test results
- Compact, ergonomic design
- Memory, stores up to 200 results

ACON works strategically with its sister companies to produce high quality medical products that meet and exceed the customers exact quality expectations.

We are committed to continuously monitoring our quality systems to ensure their effectiveness and compliance to all applicable regulatory and quality requirements. Exceeding customer expectations and maintaining customer loyalty are our top priorities.

ACON LABORATORIES, INC.

Your**Trusted**Partner



ACON Laboratories, Inc.

5850 Oberlin Drive, #340, San Diego, CA 92121, U.S.A.

Tel: 1.800.838.9502

Fax: 1.858.200.0729

Email: info@aconlabs.com