

Convenience, Accuracy
and Reliability

On·Call® Extra Mobile
Blood Glucose Monitoring System

User-Friendly Features for Easier Diabetes Management

The On Call® Extra Mobile Blood Glucose Monitoring System offers user friendly features for easier diabetes management at competitive prices. The smart and powerful On Call® Extra Mobile Meter is compatible with select remote patient monitoring platforms for real-time data management.

On Call Extra Mobile is the perfect solution to glucose monitoring when convenience matters.






On·Call® Extra Mobile

- Bluetooth enabled for wireless data management
- Now compatible with the diasend® mobile apps for iOS and Android
- Large display size for easy reading of test results
- Convenient size for on-the-go testing






diasend.
by *glooko*



Reliable and Accurate

-  Clinically Proven Accuracy in compliance with EN ISO15197:2015
-  Hypo, Hyper Alerts
-  500 test memory with date and time

Easy to Use

-  4 Second Testing
-  Only 0.4 μ L Blood Sample
-  Re-apply a second drop of blood within 3 seconds to avoid wasted strips
-  No Coding Required
-  Audible fill detection

Convenience, Accuracy and Reliability

On·Call® Extra Mobile

Blood Glucose Monitoring System

Specifications

FEATURE	SPECIFICATION
Technology	Biosensor/Electrochemical, Glucose Oxidase (GOD)
Measurement Range	0.6 – 33.3 mmol/L (10 – 600 mg/dL)
Units of Measure	Millimoles per liter (mmol/L) or milligrams per deciliter (mg/dL)
Sample Type	Fresh whole blood (Capillary)
Hematocrit Range	25 – 60%
Minimum Sample Size	0.4 µL
Test Time	4 seconds
Memory Storage	500 results with time and date
Battery	One (1) 3.0V CR 2032 Lithium coin cell battery
Battery Life	500 measurements
Meter Size (L x W x H)	75 mm x 40 mm x 13 mm
Display Size (L x W)	50 mm x 25 mm
Weight	36.0 g (with battery installed)
Operational Environment Conditions	5 – 45 °C (41 – 113 °F) and 10 – 90% relative humidity
Data Transfer	Bluetooth 4.2

Ordering Information

System and Meter

PRODUCT NAME	CATALOG NO	METER	TEST STRIPS	CONTROL SOLUTION 1	LANCETS	LANCING DEVICE*
On Call Extra Mobile Blood Glucose Monitoring System	G115-10T	•	10	•	10	•
On Call Extra Mobile Blood Glucose Meter	G115-11T	•				•
	G115-12T	•		•	10	•
	G115-13T	•		•		

All Systems and Meter Kits are available with a manual, warranty card, carrying case and Quick Reference Guide.

*Lancing device comes with clear cap for testing on forearm and palm

Strips and Accessories

PRODUCT NAME	CATALOG NO	CONTENTS
On Call Extra Blood Glucose Test Strips	G135-10H	Available in different configuration of 25 Test Strips per vial and 50 Test Strips per vial, Package Insert
	G135-11H	Available in quantities of 25 and 50 Individually foil wrapped Test Strip, Package Insert
On Call Extra Glucose Control Solution	G125-12H	1 Control Solution 0, 1 Control Solution 1, and 1 Control Solution 2; Package Insert
On Call Lancets	G124-10A	100 Lancets (25/bag)
On Call Lancing Device	G124-11A	1 Lancing Device, Package Insert



aconlabs.com

ACON Laboratories, Inc.
San Diego, CA 92121, U.S.A.
Tel: 1.800.838.9502
Fax: 1.858.200.0729
Email: info@aconlabs.com